



---

SIMONS ELECTRON  
MICROSCOPY CENTER

---

NEW YORK STRUCTURAL BIOLOGY CENTER



**Simons Electron Microscopy Center EM course**


## Challenges in Biological EM & Sample Prep

**2017/01/18**



# RT & Cryo Sample Prep Methods




adapted from  
Pilhofer et al.,  
MCB 2010





# Negative Staining

- Heavy metal salt solution surrounds sample
- Continuous carbon support film
- Protocol: glow discharge, sample, wash, stain
- SEMC: UA/UF, PTA, ammonium molybdate
- **Advantages:** high contrast, easy to learn, high SNR, radiation resistant, 3D reconstruction possible
- **Disadvantages:** structural collapse & flattening artifacts, non-native environment, ~20 Å max resolution




Baker, 2007




# Plunge Freezing

- Sample suspended in physiological buffer
- Holey carbon support film: C-flats, Quantifoil
- Protocol: glow discharge, sample, blot, plunge freeze
- SEMC: Gatan CryoPlunge Freezer 3, FEI Vitrobot, manual plunge freezer
- **Advantages:** no fixation/dehydration/staining artifacts, learning curve, random orientation, higher resolution than stain
- **Disadvantages:** low contrast, low SNR, radiation sensitive, difficult to visualize <100 kD, freezing artifacts



Protochips.com



[www.mcb.ucdavis.edu/cryoem/microscopy101.html](http://www.mcb.ucdavis.edu/cryoem/microscopy101.html)




# FIB/SEM vs Thin Section Sample Prep

- Chemical fixation
- Staining
  - En bloc, enhanced contrast and electrical conductivity
- Dehydration
- Embedding
- Au/Pd coat
  - Conductivity




- Chemical fixation
- Dehydration
- Embedding
- Sectioning
- Staining

Cryofixation: High pressure freezing  
Dehydration: Freeze substitution





# Cryo FIB Milling









## Traditional substrates for cryo-EM

- Proteins interact with surfaces present during the blotting process
  - ➡ Denaturation of proteins, preferential orientations
- Electron radiation induces motion of the particles and substrates
  - ➡ Image blurring
- Additional layer of carbon reduces signal to noise per particle
  - ➡ alignment more difficult
- Overall lack of reproducibility from grid to grid



# Traditional substrates for cryo-EM



Quantifoil, C-flat  
Cryomesh

Passmore, 2014



# Gold Grids


- Holey gold foil on gold mesh grid

## Advantages:

- Prevents differential thermal contraction when freezing
- Reduces beam-induced specimen movement
- Combined with direct detector technology allows for near atomic resolution

## Disadvantages:


- Difficult to find focus due to lack of amorphous substrate



Russo & Passmore, 2015




# Gold Grids



Russo & Passmore, 2015



## Current CryoTEM Specimen Preparation









**Spotiton**




**Typhon**



## Obtaining a 3D structure from a 2D image




From W. Baumeister et al. [Trend in Cell Biology 9\(1999\)81](#)






## Obtaining a 3D structure from a 2D image













# Vitrification process for CryoTEM




Slide courtesy Gabriel Lander
















Slide courtesy Gabriel Lander






Slide courtesy Gabriel Lander






Slide courtesy Gabriel Lander






Slide courtesy Gabriel Lander






Slide courtesy Gabriel Lander





Slide courtesy Gabriel Lander





Slide courtesy Gabriel Lander